

# ECOKING

Soft and noiseless non-PVC film

Ideal for:



Urine and drainage applications

STERILE R

WIPAK is a member of Sterile Barrier Association



# ECOKING

Soft and noiseless non-PVC film



## Material Composition

- Thickness 80, 150, 180 microns

## Applications

- For urine and drainage bags
- For irradiation sterilization
- Heat sealable

## Benefits

- Soft material
- Noiseless
- No PVC
- Available also in white color

Wipak reserves right to make changes at any time without notice.

WIPAK Oy  
P.O. Box 45  
Wipaktie 2  
FI-15561 Nastola, Finland  
Tel: +358-20510 311  
Fax: +358-20510 3333  
email: [ind@wipak.com](mailto:ind@wipak.com)  
[www.wipak.com](http://www.wipak.com)

The WIPAK logo, consisting of a stylized green 'W' followed by the word 'WIPAK' in a bold, green, sans-serif font.